

Manual BerryKing Healthbulb

Technical specifications:

- Use: 6 hours
- Implementation of standards: GB4706.48-2000
- Product type: humidifier
- Water capacity: 400 ml
- Material: ABS / PP / PET
- Power: 2W
- Size: 153 x 88 x 88 mm
- Input voltage: DC 5V
- Cable length: 120 cm (USB interface)

Features:

- Press briefly: Switch on / off
- Press for 2 seconds: Switch on the light
- Press for 5 seconds: View light



Application:

- ① : Open the cover and separate it from the sleeve
- ② : Insert the spring
- ③ : Insert cotton swabs (moisten before use for the first time)
- ④ : Reattach the lid with the rod cover
- ⑤ : Place decoration in the lamp (if desired) and fill with water (max. 400ml)
- ⑥ ⑦ : Close the lamp with the lid including cotton swabs
- ⑧ : Connect the USB cable to the PC or the socket (adapter required) and switch it on



In the box:

- BerryKing Healthbulb
- Decoration (stones, palm trees)
- Cotton swab
- Manual